

## STENTS (Bare metal and drug-eluting)

### Cobalt chromium alloy ASTM F75

Metal	%
Chromium	28.5
Molybdenum	6
Nickel	0.25
Manganese	0.5
Tungsten	<0.2
Cobalt	Balance

### Stainless steel 316L

Metal	% average
Manganese	2
Chromium	17
Nickel	12
Molybdenum	2.5
Iron (not tested in MELISA)	Balance

### Nitinol

Metal	% average
Nickel	50
Titanium (dioxide and sulphate)	50

### Platinum chromium

Metal	% average
Platinum	33
Chromium	18
Molybdenum	3
Manganese	1
Iron (not tested in MELISA)	Balance